




## Sample rack, 28 x 50 mL, (PVC)

6.2041.490

Sample rack for 919 Autosampler plus, 858 Professional Sample Processor suitable for 28 sample vessels. The following sample vessels can be used: 6.2747.010.

Plastic: Polyvinyl chloride (PVC)

### Phụ kiện tùy chọn

| Order no.  | Description         |
|------------|---------------------|
| 6.2747.010 | Sample vessel 50 mL |

With blue screw cap and mL graduation. Set of 25 pieces.



6.1458.080

Transfer head for 919 IC Autosampler plus for VA and CVS applications.

