




Sample rack, 36 x 6 mL, (Al)

6.2041.720

Sample rack for 874 Oven Sample Processor, suitable for 36 sample vessels. The following sample vessels can be used: 6.2419.000.

Material: Aluminum (Al)

Phụ kiện tùy chọn

Order no.	Description
6.2419.000	Sample glass, 6 mL, 1000 pieces
	Sample glass for up to 6 mL sample, suitable for 6.1448.050 septum, 1,000 pieces



6.1448.050

Aluminum septum caps / 1000 pieces

Including Al sealing.

