



Sample rack insert, 54 x 75 mL, for water bath, (Al)

6.2041.900

Sample rack insert for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 54 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400.

Material: Aluminum (Al)

Phụ kiện tùy chọn

Order no.	Description
6.1432.210	Sample beaker, clear glass, 75 mL, 1 piece
	Sample beaker, clear glass, 75 mL, suitable for sample racks: 6.2041.340, 6.2041.350, 6.2041.800 and 6.2041.900, 1 piece



6.1462.050 Robotic arm with titration head for 786 Swing Head, left or right swinging
Titration head for Swing Head



6.1462.150 Robotic arm with titration head and beaker sensor for 786 Swing Head,
left or right-swinging
Titration head with beaker sensor for 786 Swing Head, 2 x M10, 1 x SGJ
14.5, 2 x buret tip



6.1459.400 Sample beaker, plastic (PE), 75 mL, 144 pieces

Sample beaker, plastic (polyethylene, PE) 75 mL, suitable for sample
racks: 6.2041.340, 6.2041.350, and 6.2041.800, suitable for 6.02041.040
OMNIS Software sample rack, 144 pieces.



6.2840.000 Water bath for Robotic Sample Processors and Robotic Titrosampler
Water bath rack for Robotic Sample Processors and Titrosampler.



6.2711.070

Drip pan for 789/815/855

Drip pan for 789, 815, 855 Robotic Sample Processors XL

