




Bottle holder to Extension Module - Liquid Handling

6.2057.200

Phạm vi giao hàng 6.2057.200

Qt.	Order no.	Description
1 PCS	6.2621.030	Hexagon key 4 mm



Phụ kiện tùy chọn

Order no.	Description	
6.1562.160	Transfer tubing 15 mL, 2 x M6 For liquid handling in ion chromatography. To avoid contamination of the Dosino.	
6.1608.100	PE bottle 50 mL Bottle for Extension Module Liquid Handling.	
6.2762.000	Mixing vessel 50 mL to Extension Module - Liquid Handling	