



Scale pan for sample rack 6.2068.010

6.2067.100

The scale pan is designed for the 120 mL PP beaker (6.1459.300) and is placed into the sample rack of the 864 Balance Sample Processor.

Phụ kiện tùy chọn

Order no.	Description
6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces
	Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces

