




Sample rack, 20 x 120 mL, (PP)

6.2068.010

Sample rack for 864 Balance Sample Processor, suitable for 20 sample beakers. The following sample beakers can be used: 6.1459.300.

Plastic: Polypropylene (PP)

Phụ kiện tùy chọn

Order no.	Description
6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces
	Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces



6.1562.240

Pipetting tips 1-10 mL (200 pieces)

Eppendorf tips for use with 864 Robotic Balance Sample Processor



6.2067.100

Scale pan for sample rack 6.2068.010

The scale pan is designed for the 120 mL PP beaker (6.1459.300) and is placed into the sample rack of the 864 Balance Sample Processor.

