



Sample rack, 24 x 120 mL, (PVC)

6.2068.020

Sample rack for 815 Robotic Soliprep, suitable for 24 sample beakers. The following sample beakers can be used: 6.1459.300, 6.1459.307.

Plastic: Polyvinyl chloride (PVC)

Phụ kiện tùy chọn


Order no.	Description
6.1459.300	Sample beaker, plastic (PP), 120 mL, 250 pieces
	Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 250 pieces



6.2067.000 Insert for 30x11 mL for Soliprep sample rack
Sample rack insert for 30x11 mL (6.2743.050)



6.2067.010 Insert for 24xSyringe needles for Soliprep sample rack
Sample rack insert for 24xsyringe needles (luer connection)



6.2067.020 Insert for 12x30 mm syringe filters for Soliprep Sample rack
Sample rack insert for 12x30 mm syringe filters (luer connection)



6.1459.307 Sample beaker, plastic (PP), 120 mL, 1 piece

Sample beaker, plastic (polypropylene, PP), 120 mL, suitable for sample racks: 6.2041.470, 6.2041.840, 6.2041.850, 6.2068.010, 6.2068.020, 6.2068.030, and 6.9920.141, 1 piece

