



Power supply unit for 949 pH Meter

6.2166.200

Tabletop power supply unit, 100–240 V, 50–60 Hz, 5 V DC, for 949 pH Meter

Efficiency Level VI

Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

Discover the 949 pH Tabletop Meter – Your compact pH measuring instrument for quick and reliable laboratory measurements.

- Complete pH measuring instrument with color display
- Precise pH, mV and ORP measurements with integrated temperature display
- Simple calibration with Metrohm, USA, and NIST buffers up to 3 points or by 2-point calibration with in-house buffers
- Clearly structured display of the calibrated values using symbols for simple read-off
- 3 stability criteria for measured value recognition for reliable results
- GLP conformity of the measured values with embedded date and time
- Overview of calibration data with monitoring function
- RS232 port enables connection of a printer.



Included in the package are the 949 pH Meter, the electrode holder, and the power supply unit for immediate use in the laboratory.
