



# Buffer solution pH 7.00 (30 x 30 mL)

6.2307.210

Buffer solution pH 7.00 (25 °C) in single use sachets, colorless

Phụ kiện tùy chọn

| Order no.  | Description                                                                                                                                 |
|------------|---------------------------------------------------------------------------------------------------------------------------------------------|
| 6.2307.230 | Buffer solutions pH 4, 7 and 9<br>Mixed buffer solutions pH 4.00/7.00/9.00 (25 °C) in single use sachets, colourless, box of 3 x 10 x 30 mL |



6.2307.200 Buffer solution pH 4.00 (30 x 30 mL)

Buffer solution pH 4.00 (25 °C) in single use sachets, colorless



---

6.2307.220 Buffer solution pH 9.00 (30 x 30 mL)

Buffer solution pH 9.00 (25 °C) in single use sachets, colorless



---

6.2307.100 Buffer solution pH 4 (500 mL)

Ready-to-use buffer solution pH = 4 (25 °C), colored, with closure seal



---

6.2307.110 Buffer solution pH 7 (500 mL)

Ready-to-use buffer solution pH = 7 (25 °C), colored, with closure seal



6.2307.120

Buffer solution pH 9 (500 mL)

Ready-to-use buffer solution pH = 9 (25 °C), colorless, with closure seal



---

6.2325.000

pHit kit

Maintenance kit for pH electrodes

The kit contains:

- 50 mL cleaning solution
- 50 mL 3M KCl solution
- 50 mL storage solution
- 2 Storage vessels
- Instructions for use

