



Capillary cutter

6.2621.080

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.

Phụ kiện tùy chọn

Order no.	Description
6.1803.040	PTFE capillary 0.5 mm i.d. / 1 m Capillary for sample handling in IC.



6.1803.080

PTFE capillary 0.97 mm i.d. / 1 m

Capillary for filtration, for IC Filtration Sample Processor, IC Liquid Handling Filtration Unit



6.1831.010

PEEK capillary 0.25 mm i.d. / 3 m

For all IC components.



6.1831.030

PEEK capillary 0.75 mm i.d. / 3 m

For all IC components.



6.1831.060

PEEK capillary 0.5 mm i.d., 1 m

For inline dialysis and for Sample Changers. For IC Dialysis Unit, IC Sample Processor, Compact IC Autosampler, IC Filtration Sample Processor, IC Dialysis Sample Processor, IC Liquid Handling Dialysis Unit

