



## Filtration membrane reg. cellulose, 0.45 µm, 50 pieces

6.2714.050

Nominal pore size: 0.45 µm. Set of 50 pieces. For inline filtration.


### Phụ kiện tùy chọn

| Order no.   | Description                                                                                                                                                                                               |
|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 6.02729.110 | Ultrafiltration cell 2                                                                                                                                                                                    |
|             | For the Inline ultrafiltration of samples together with the 858 Professional Sample Processor / 919 IC Autosampler plus. As an alternative to the time-consuming manual filtration with a syringe filter. |



6.05330.010 IC equipment: Inline Ultrafiltration 2 - pull mode

Accessory set for Inline Ultrafiltration 2 - pull mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.




---

6.05330.110 IC equipment: Inline Ultrafiltration 2 - push mode

Accessory set for Inline Ultrafiltration 2 - push mode. For use with the 858 Professional Sample Processor/ 919 IC Autosampler plus.

Recommended set for the combination with Inline Dilution (MIDT).



---

6.05330.210 IC equipment: Inline Ultrafiltration 2 - MiPT

Accessory set to combine Inline Ultrafiltration 2 and MiPT. For use with the 858 Professional Sample Processor.



---

6.2831.010 Plastic tweezers

For dialysis and filtration membranes

