




Dialysis cell

6.2729.100

With 4 x 6.2744.060 Blind stoppers. For inline dialysis

Phạm vi giao hàng 6.2729.100

Qt.	Order no.	Description
4 PCS	6.2744.060	Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Phụ kiện tùy chọn

Order no.	Description
6.2744.000	Pressure screw PVDF With UNF 10/32 connection. For IC instruments, inline dialysis (connection of PTFE capillaries)

