



Septum stopper

6.2730.020

With E.301.0041 O-ring. For 6.1448.010 Septum and 6.1414.030 Titration vessel lid. Used for KF applications.

Phạm vi giao hàng 6.2730.020

Qt.	Order no.	Description
1 PCS	6.2730.010	Screw nipple
		Made of PP. For 6.1448.010 Septum or for mounting 6.1403.040 Drying tube or 6.0338.100 Electrode in 6.1414.030 and 6.1414.040 KF Titration vessel lids. Used for KF titrations.



1 PCS

6.1448.010

Septum 12 mm, 5 pieces

Set of 5 items.

