




EMPTY sensor for level sensor for 2 L bottles

6.2769.110

This sensor, mounted in the 6.1626.000 bottle attachment, results in the level sensor EMPTY for 2 L bottles. Included are two sensor rods, one aspiration tubing and two holders for fixing the tubing to one of the sensor rods.

Phạm vi giao hàng 6.2769.110

Qt.	Order no.	Description
1 PCS	6.1819.100	FEP aspiration tubing to 2 L bottle

Aspiration tubing to bottle attachment with M6 connector. For the production of eluents.

