



Aspiration filter

6.2821.090

Pore size: 27 μm . Set of 5 pieces. For aspiration tubing (6.1834.000) and filter tubes (6.1821.040) and (6.1821.050).

Phụ kiện tùy chọn

Order no.	Description
6.02744.050	Adapter Luer inner / M6 outer

Adapter for syringe for pre-rinsing the aspiration filter (6.2821.090).
The adapter can be screwed onto the aspiration filter. A syringe can be used to rinse the aspiration filter with ultrapure water or eluent.



6.2816.020

Syringe 10 mL with Luer connection

For various applications in IC and VA

