



Standard pump head PEEK

6.2824.110

Standard pump head for intelligent IC instruments

Phạm vi giao hàng 6.2824.110

Qt.	Order no.	Description
2 PCS	6.2741.020	Piston seal (orange) For all standard pump heads



1 PCS

6.2824.170

Inlet valve PEEK

Metal-free check valve for the inlet side of the IC pump heads of the devices 940, 930, 925, 883, 882, 881, and 850.



1 PCS

6.2824.160

Outlet valve PEEK

Metal-free check valve for the outlet side of the IC pump heads of the devices 940, 930, 925, 883, 882, 881, and 850.

