



Water bath for Robotic Sample Processors and Robotic Titrosampler

6.2840.000

Water bath rack for Robotic Sample Processors and Titrosampler.

Phạm vi giao hàng 6.2840.000

Qt.	Order no.	Description
2 PCS	6.1801.130	PVC tubing / 8 mm / 12 mm / 2 m



Phụ kiện tùy chọn

Order no.	Description	
6.2041.900	Sample rack insert, 54 x 75 mL, for water bath, (Al) Sample rack insert for 815 Robotic USB Sample Processor XL, 855 Robotic Titrosampler "Basic", suitable for 54 sample beakers. The following sample beakers can be used: 6.1432.210, 6.1459.400. Material: Aluminum (Al)	
6.2058.030	Head extension 10 mm for Sample Processors Head extension for use with the waterbath rack.	