



Rinsing station for pipetting equipment

6.2841.110

Rinsing station for pipetting equipment

Phạm vi giao hàng 6.2841.110

Qt.	Order no.	Description
2 PCS	6.1805.120	FEP tubing / M6 / 100 cm
		With light and kink protection



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



2 PCS

6.1808.090

Thread adapter M8 outer / M6 inner

Outer thread M8, inner thread M6, e.g. for mounting 6.2418.XXX Air tubes in 6.1429.XXX Reaction vessels of stability measuring instruments.



1 PCS

6.1816.080

Silicone tubing 6 mm i.d. / 2 m

Drainage tubing to the "Waste Collector".

