




# HT reactor plate, complete, for Professional Reactor

6.2845.110

Reactor plate with reactor for PCR applications at high temperatures ( $> 120\text{ }^{\circ}\text{C}$ ). With PEEK reactor capillaries.


Phạm vi giao hàng 6.2845.110

| Qt.   | Order no.  | Description      |
|-------|------------|------------------|
| 3 PCS | 6.2744.060 | Threaded stopper |

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



Phụ kiện tùy chọn

| Order no.  | Description                                              |                                                                                     |
|------------|----------------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2057.240 | Holder for reactor 62845xxx                              |                                                                                     |
|            | Holder for fastening a reactor in the column thermostat. |  |

---