



Exchange Unit 1 mL

6.3026.110

Exchange Unit with integrated data chip with 1 mL glass cylinder and light protection. PCTFE /PTFE flat stopcock, FEP tubing connection, antidiffusion buret tip and standard amber glass reagent bottle.

Phạm vi giao hàng 6.3026.110

Qt.	Order no.	Description
1 PCS	6.1543.060	Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



1 PCS

6.1805.170

PTFE capillary / M6 / 18 cm


With light protection and protection against kinks.



Phụ kiện tùy chọn

Order no.	Description
-----------	-------------

6.1576.110 Exchange Unit base plate, 1 mL



6.1608.030 Round glass bottle / 1000 mL / GL 45



6.1608.040 Bottle PE / 1000 mL / GL 45
For Exchange Units. Bottle for auxiliary solutions.

