



Exchange Unit 20 mL

6.3026.220

Exchange unit with integrated data chip with 20 mL glass cylinder and light protection. PCTFE /PTFE flat stopcock, FEP tubing connection, anti-diffusion buret tip and standard amber glass reagent bottle.

Phạm vi giao hàng 6.3026.220

Qt.	Order no.	Description
1 PCS	6.1543.060	Tip / M6 thread

Tip with M6 thread. Together with the 6.1805.160 capillary and the 6.1446.030 link stopper, it makes up the entire 6.1537.010 buret tip.



Phụ kiện tùy chọn

Order no.	Description	
6.1563.040	Thermostat jacket for 806 Exchange Units For 6.3026.XXX Exchange Units	
6.1608.030	Round glass bottle / 1000 mL / GL 45	
6.1608.040	Bottle PE / 1000 mL / GL 45 For Exchange Units. Bottle for auxiliary solutions.	
6.1576.220	Exchange unit base plate / 20 mL	