



Dosing unit PCR 5 mL

6.3033.150

5 mL dosing unit specially suitable for conveying solutions for post-column derivatization (PCR) in IC.

Phạm vi giao hàng 6.3033.150

Qt.	Order no.	Description
1 PCS	6.2739.000	Wrench For tightening connectors



1 PCS

6.2803.010

Grease (2 g)

Special quality (without silicone). For grease-free ground joint connections see 6.2713.XXX



Phụ kiện tùy chọn

Order no.	Description	
6.2057.250	Dosino holder for PCR	
	Dosino holder for IC instruments for PCR Dosino.	

