



Accessory kit: Vario/Flex Basic

6.5000.000

Basic set of accessories for a Professional IC Vario or Compact IC Flex.

Phạm vi giao hàng 6.5000.000

Qt.	Order no.	Description
1 PCS	6.01816.010	Silicone tubing, 8/12 mm, 2 m
		Drainage tubing



1 PCS

6.1602.210

Attachment 5 x UNF 10/32 for «Waste Collector»

Attachment for holding up to 5 waste capillaries.



1 PCS

6.1807.010

Y-connector for tubing 6-9 mm i.d.

Connector for waste tubing



1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.



2 PCS

6.1816.020

Silicone tubing 6 mm i.d., 1 m

For waste tubing.



1 PCS

6.2057.220

Holder to Waste Collector

Holder for fastening the Waste Collector to the IC instrument.



1 PCS

6.2151.020

Cable USB A - USB B 1.8 m

USB connecting cable



1 PCS

6.2251.000

Colored sleeves for capillaries

Colored pieces of heat shrink tubing for capillary coding. Three pieces each of five different colors.



1 PCS

6.2617.010

Tool for piston seal

For removing and assembling the piston seal for all standard pump heads



2 PCS 6.2621.000 Adjustable wrench
Max. opening: 20 mm. For IC instruments



1 PCS 6.2621.030 Hexagon key 4 mm



1 PCS 6.2621.050 1/4 in. wrench
For 1/4 in. screws. For IC instruments



1 PCS 6.2621.080 Capillary cutter

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8", including 5 replacement blades.



1 PCS

6.2621.100

3 mm hex key for IC Sample Processors



2 PCS

6.2739.000

Wrench

For tightening connectors



1 PCS

6.2743.080

Stopper for overflow (ProfIC), 5 pieces

For Professional IC instruments



1 PCS

6.2765.000

Container with GL 45 drainage to the "Waste Collector"

Container with drainage to the "Waste Collector".

