




Differential amplifier plug F / 230 V EU

6.5104.030

Differential amplifier with power supply.

Phạm vi giao hàng 6.5104.030

Qt.	Order no.	Description
1 PCS	6.2117.000	Power adapter EU 230 V / 9 V AC
With EU mains plug. Input: 220...240 V AC. Output: 9V AC/0.2 A. Used with 587 Conductometer, 588 pH Meter, 627 O ₂ Meter		



Differential amplifier with plug F, for Titrados (power supply unit must be ordered separately: 6.2117.000 power unit EU 230 V / 9 V AC)

Also available as a complete set: 6.5104.030 differential amplifier / plug F / 230 V EU

