



IQ/OQ Kit for IC

6.5333.000

The IQ/OQ Kit contains all parts and standard solutions required for IQ/OQ in ion chromatography.

Phạm vi giao hàng 6.5333.000

Qt.	Order no.	Description
1 PCS	6.1803.050	PTFE capillary 0.5 mm i.d. / 20 cm Capillary for inline dialysis and filtration, for Dialysis Unit, IC Filtration Sample Processor



1 PCS

6.1831.170

OQ capillary

PEEK capillary used in OQ procedures for Professional IC equipment



1 PCS

6.2322.030

OQ Standard Solutions for IC

Part of the IQ/OQ Kit for IC (6.5333.000). Cannot be ordered individually.



1 PCS

6.2744.040

2 x UNF 10/32 coupling

For connecting 1/16 in. capillaries. For IC instruments



1 PCS

6.2744.060

Threaded stopper

For UNF 10/32. Stopper for IC, e.g. for the sealing of columns.



1 PCS

6.2744.070

Pressure screw short

Short version. With UNF 10/32 connection. 5 pieces. For the connection of PEEK capillaries

