



Equipment with 2x 807 Dosing Units for VA and CVS measurements

6.5339.500

Accessory set for 2 auxiliary solutions for expanding VA or CVS systems. Contains 2x 807 Dosing Units 2 mL, 2x 100 mL clear glass bottles, bottle holder for 2 bottles, and connector material. Without Dosinos.

Phạm vi giao hàng 6.5339.500

Qt.	Order no.	Description
2 PCS	6.1608.050	Drying bottle / 100 mL / GL 45



2 PCS

6.1618.060

Thread adapter GL 45 to GL 45

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



2 PCS

6.1805.530

FEP tubing / M6 / 2 m

With light and kink protection



2 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



1 PCS

6.2055.110

Bottle holder or Dosino holder for two bottles
/Dosinos

Holder for two 250 mL bottles or two Dosinos



807 Dosing Unit with integrated data chip with 2 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Phụ kiện tùy chọn

Order no.	Description
2.800.0010	800 Dosino
	The 800 Dosino is a drive with write/read hardware for intelligent Dosing Units. With fixed cable (length 150 cm).

