



Basic equipment with 2x 807 Dosing Units for VA measurements for 884 Professional VA

6.5339.510

Accessory set for 2 auxiliary solutions for expanding manual VA systems. Includes one 807 Dosing Unit 2 mL and one 807 Dosing Unit 5 mL, clear glass bottles 100 mL, 250 mL and 2 L, bottle holders for 2 bottles, one 4-way micro tip, and connector material. Without Dosinos.

Phạm vi giao hàng 6.5339.510

Qt.	Order no.	Description
4 PCS	6.1446.040	Threaded stopper / M6



1 PCS

6.1602.190

Eluent Bottle Head, 1 x SGJ 14, 4 x M6

Eluent bottle head to connect up to 4 tubings with M6 connectors. For eluent bottle of the Eluent synthesizer.



1 PCS

6.1608.050

Drying bottle / 100 mL / GL 45



1 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS

6.1608.090

Bottle / 250 mL / GL 45



2 PCS

6.1618.060

Thread adapter GL 45 to GL 45

Thread adapter GL45 to GL45; e.g. for reagent bottles with extended flange



1 PCS

6.1805.120

FEP tubing / M6 / 100 cm

With light and kink protection



2 PCS

6.1805.530

FEP tubing / M6 / 2 m

With light and kink protection



4 PCS

6.1808.000

Connection bushing M6



1 PCS

6.1815.010

Spiral band 0.5 m

For holding together different cables or tubing.



2 PCS

6.1819.010

PTFE tube

As gas inlet in washing bottle of VA Stand



1 PCS

6.1824.000

Four-way micro dosing tip

With 4 x M6 thread. For VA measuring stands.



1 PCS

6.1829.070

Tubing set for Professional VA/CVS instruments

Tubing set for the measuring head of Professional VA/CVS instruments. Comprised of 4 FEP tubings with M6 connector, tapering tip and protected against kink (2 x 82 mm long und 2 x 110 mm long).



1 PCS

6.2055.110

Bottle holder or Dosino holder for two bottles
/Dosinos

Holder for two 250 mL bottles or two Dosinos



1 PCS

6.2621.080

Capillary cutter

Cutting tool for simple cutting of plastic capillaries, 1/16" and 1/8",
including 5 replacement blades.



1 PCS

6.3032.120

807 Dosing Unit 2 mL

807 Dosing Unit with integrated data chip with 2 mL glass cylinder and
light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass
thread. FEP tubing connection, antidiffusion tip.



1 PCS

6.3032.150

807 Dosing Unit 5 mL

807 Dosing Unit with integrated data chip with 5 mL glass cylinder and
light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass
thread. FEP tubing connection, antidiffusion tip.

