



Equipment: 807 Dosing Unit for pipetting system for Professional VA/CVS instruments

6.5339.620

Accessory kit for fully automated CVS pipetting system in combination with the 858 Professional Sample Processor. Contains one 807 Dosing Unit 5 mL and 50 mL each, clear glass bottles 2 L, pipetting tubing, rinsing station, safety shield and other accessories. Without Dosinos.

Phạm vi giao hàng 6.5339.620

Qt.	Order no.	Description
1 PCS	6.1456.210	Measuring vessel / 10-90 mL
		Measuring vessel for automated VA systems. Dilution vessel for Advanced IC Sample Processor



1 PCS

6.1462.270

Robotic arm with transfer head, for 786 Swing Head,
right-swingingRobotic arm for pipetting system with 858 Professional Sample Processor
VA and CVS

1 PCS

6.1562.120

Pipetting tubing / 10 mL with holder for sample
processors

Pipetting tubing for liquid transfer.



1 PCS

6.1602.115

Bottle attachment / S 40 / green

For direct dosing from reagent bottles with S 40 thread (Merck...).



2 PCS

6.1608.070

Eluent bottle / 2 L / GL 45



1 PCS 6.1805.080 FEP tubing / M6 / 25 cm
With light and kink protection.



1 PCS 6.1805.120 FEP tubing / M6 / 100 cm
With light and kink protection



3 PCS 6.1805.530 FEP tubing / M6 / 2 m
With light and kink protection



1 PCS 6.1808.010 Coupling piece 2 X M6
2 M6 outer threads. Part of the screw connections for Exchange Units and
676 Sample Changer



1 PCS

6.1808.050

Tubing adapter olive / M8 outer

1 M8 outer thread and 1 tubing olive. E.g. for thermostat jacket of exchange units and stability measuring instruments.



1 PCS

6.1808.060

T connector / M6 inner

3 M6 inner threads.



1 PCS

6.1819.080

Cannula for gas feeding

For stability measuring instruments.



1 PCS

6.1823.010

Guide shaft to pipetting tubing

For assembly in the transfer robotic arm on the Robotic Sample Processors



2 PCS

6.1829.080

Aspiration tubing / M6 / 0.3 m




1 PCS

6.2622.040

Positioning plate

Positioning plate for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.2622.030

Safety shield accessories

Accessories for pipetting system with 858 Professional Sample Processor VA and CVS. Contains stand and fastening material for the 6.2751.200 safety shield and the 6.2622.020 positioning plate.



2 PCS

6.2709.070

Guiding sleeve

For fixation of buret and dosing tips in titration heads and electrode holders



2 PCS

6.2744.200

Adapter UNF 10/32 outer / M6 inner

Connecting 1/16 in. PEEK capillaries and transfer tubing (M6)



1 PCS

6.2751.200

Safety shield

Safety shield for pipetting system with 858 Professional Sample Processor VA and CVS



1 PCS

6.2841.110

Rinsing station for pipetting equipment

Rinsing station for pipetting equipment



1 PCS

6.3032.150

807 Dosing Unit 5 mL

807 Dosing Unit with integrated data chip with 5 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Phụ kiện tùy chọn

Order no.	Description
2.800.0020	800 Dosino Drive with write/read hardware for intelligent dosing units. With permanently attached cable (length 0.65 m).