



Equipment for determination of temperature correction with Rancimats and PVC Thermomats.

6.5616.100

Set for exact temperature adjustment

Phạm vi giao hàng 6.5616.100

Qt.	Order no.	Description
1 PCS	6.1111.010	Pt100 temperature sensor for stability measuring devices (17.5 cm)

Pt100 temperature sensor (length 17.5 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with a 711 Liquino and instruments for stability measurement (743 Rancimat, 763 PVC Thermomat, 892 Professional Rancimat, or 895 Professional PVC Thermomat).



2 PCS

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments



1 PCS

6.2326.000

Silicone oil for stability measuring instruments (50 mL)

Silicone oil for the determination of temperature correction in stability measuring instruments



1 PCS

6.2418.000

Air tube for oil measurements

Set of 12 pieces. In combination with 6.1429.040.

