



Equipment for determination of temperature correction with Biodiesel Rancimats.

6.5616.110

Set for exact temperature adjustment

Phạm vi giao hàng 6.5616.110

Qt.	Order no.	Description
1 PCS	6.1111.020	Pt100 temperature sensor for stability measuring devices (30 cm)

Pt100 temperature sensor (length 30.0 cm) with fixed cable and mini-DIN-plug. This glass-free temperature sensor is well suited for use with instruments for stability measurement (873 Biodiesel Rancimat and 893 Professional Biodiesel Rancimat).



2 PCS

6.2042.040

Clamp for temperature sensor

Holding clamp for Pt100 temperature sensor for stability measuring instruments



1 PCS

6.2326.000

Silicone oil for stability measuring instruments (50 mL)

Silicone oil for the determination of temperature correction in stability measuring instruments

