



Equipment for reagent replacement

6.5617.000

Equipment for automatic reagent replacement with KF coulometer and connected Dosino.

Phạm vi giao hàng 6.5617.000

Qt.	Order no.	Description
1 PCS	6.1446.060	Stopper / B-14/15 / M10



1 PCS

6.1446.080

Threaded stopper / M8



1 PCS

6.1543.070

Aspiration tip / M6 thread

For the 661 Pump Unit



1 PCS

6.1602.105

Bottle attachment / GL 45 / green

For direct dosing from reagent bottles with GL 45 thread.



1 PCS

6.1608.030

Round glass bottle / 1000 mL / GL 45



1 PCS

6.1609.000

Adsorbing tube, large and bent

For filling with adsorber material.



1 PCS

6.1805.060

FEP tubing / M6 / 60 cm

With light and kink protection



1 PCS

6.1805.080

FEP tubing / M6 / 25 cm

With light and kink protection.



1 PCS

6.2023.020

Clip for SGJ 14/15

Clip for SGJ 14/16



1 PCS

6.2055.100

Bottle holder, double

Bottle holder for 2 x 1 L bottles. For 6.5617.000 aspiration equipment



1 PCS

6.2730.030

Stopper

Complete with nipple and O-ring. For 6.1414.030 titration vessel lids. For KF instruments, VA measurement stands, VA Computrace and IC Dilution.



1 PCS

6.3032.250

807 Dosing Unit 50 mL

807 Dosing Unit with integrated data chip with 50 mL glass cylinder and light protection, mountable to a reagent bottle with ISO/DIN GL 45 glass thread. FEP tubing connection, antidiffusion tip.



Phụ kiện tùy chọn

Order no.	Description	
2.800.0010	800 Dosino	
6.2811.000	Molecular sieve	
6.2811.010	Molecular sieve (US)	