



Transfer equipment, 0.5 - 2 mL, for 789/815

6.5619.000

Complete equipment for the transfer of solutions from 0.5 to 2.0 mL on the Robotic Sample Processor.

Phạm vi giao hàng 6.5619.000

Qt.	Order no.	Description
1 PCS	6.1446.180	Stopper, PTFE, SGJ 29, for transfer set
		Stopper for titration and volumetric KF titration in the external titration vessel



1 PCS

6.1448.000

Septum stopper B-14/15



2 PCS

6.1448.020

Septum 16 mm, 5 pieces

Set of 5 pcs.



1 PCS

6.1462.090

Robotic arm with Luer-Lock adapter for 786 Swing
Head, right swinging

Transfer robotic arm for liquid handling



1 PCS

6.1464.320

KF titration vessel / 80-250 mL / coulometric



1 PCS

6.1562.110

Pipetting tubing / 2 mL

Pipetting tubing for liquid transfer



1 PCS

6.1608.090

Bottle / 250 mL / GL 45




1 PCS

6.2001.100

Base plate, transfer set

Base plate for transfer set for liquids to be used with Robotic Sample Processors XL



1 PCS

6.2013.010

Clamping ring

For support rods with a diameter of 10 mm.



1 PCS

6.2701.040

GL 18 screw cap with hole

Made of PBT. With a hole for the 6.1448.020 septum. For the 6.5405.000 cell without diaphragm and the 6.1464.32X and 6.1465.320 titration vessels.



1 PCS

6.2751.120

Safety shield Transfer set to 789/815

Safety shield for Sample Processors with 786 Swing Head and transfer set.



1 PCS

6.2816.100

Needle with vent

With luer lock, for transfer equipment 6.5619.000.



1 PCS

6.2055.110

Bottle holder or Dosino holder for two bottles /Dosinos

Holder for two 250 mL bottles or two Dosinos

