




Consumable Kit Titration

6.5702.010

FEP tubing connection M6 (13, 25 and 40 cm) and antidiffusion tip for 6.3026.xxx exchange units and 6.3032.xxx dosing units.

Phạm vi giao hàng 6.5702.010

Qt.	Order no.	Description
1 PCS	6.1543.200	Titration tip / M6 thread
		Titration tip with antidiffusion valve.



2 PCS 6.1805.010 FEP tubing / M6 / 13 cm
With light and kink protection.



2 PCS 6.1805.080 FEP tubing / M6 / 25 cm
With light and kink protection.



2 PCS 6.1805.100 FEP tubing / M6 / 40 cm
With light and kink protection



1 PCS 6.1829.010 FEP aspiration tube M6, 0.25 m
With M 6 thread. For Dosing Unit



4 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200

