



Consumable Kit CVS

6.5703.020

Assembly of important expendable items for CVS applications with Professional CVS instruments.

Phạm vi giao hàng 6.5703.020

Qt.	Order no.	Description
1 PCS	6.0728.130	LL-Ag/AgCl reference electrode for CVS

Inner reference system with electrolyte $c(\text{KCl}) = 3 \text{ mol/L}$. For use with electrolyte vessel made of plastic 6.1245.010.



1 PCS

6.1204.610

2 mm platinum electrode tip for CVS

Metal electrode tip made of platinum (Pt), diameter of electrode disk 2 mm ± 0.05 mm, polished surface, shaft made of glass, M3 thread connection. Concentricity error: 0.3 mm. Temperature range: 0 - 50 °C.




1 PCS

6.1244.050

Drive belts for rotating disk electrodes

For Professional VA/CSV instruments (3 pcs.).



1 PCS

6.1824.000

Four-way micro dosing tip

With 4 x M6 thread. For VA measuring stands.

