



Glass vial, 1.5 mL, 100 pieces

6.7303.100

For use with the following sample changers: Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ) and 889 Sample Center.


Phụ kiện tùy chọn

Order no.	Description
6.2041.200	Sample rack, 48 x 11 mL, (Al) Sample rack for 889 IC Sample Center, suitable for 48 sample vessels. The following sample vessels can be used: 6.2743.100, 6.2743.107, 6.2743.110, 6.2743.117. Material: Aluminum (Al)



6.7303.110 Septum for 67303100, 100 pieces

For use with the 67303100 Glass vial on the following sample changers:
Liquid Autosampler CIC (TEI), MMS 5000 Autosampler (AJ), and 889
Sample Center.



6.7303.310 Sample rack 112 positions for liquid samples

For liquid samples.

