




Titration head with reinforced buret tips

6.9914.139

Titration head with 3 reinforced buret tips, one M10 screw connection and one opening for an assembly (Polytron PT 1300 D).

Phạm vi giao hàng 6.9914.139

Qt.	Order no.	Description
3 PCS	6.1808.000	Connection bushing M6



3 PCS

6.2726.090

Antidiffusion valve

For antidiffusion buret tips 6.1541.XXX, 6.1543.130, 6.1543.140,
6.1543.150, 6.1543.160, 6.1543.200

