



Sample rack, 35 x 60 mL and 3 x 250 mL (PVC)

6.9920.191

Sample rack for 858 Professional Sample Processor, 815 Robotic USB Sample Processor XL, suitable for 35 sample vessels and 3 rinsing bottles. The following sample vessels and rinsing bottles can be used: 6.1608.100, 60 mL and 250 mL Nalgene FEP vessels.

Plastic: Polyvinyl chloride (PVC)

Phụ kiện tùy chọn

Order no.	Description
6.1608.100	PE bottle 50 mL
	Bottle for Extension Module Liquid Handling.

