



# 用于 864 Robotic Balance Sample Processor 的离子交换器， 115V

2.136.0305

用于 864 Robotic Balance Sample Processor 的离子交换器（HAUG），带 115V 供电单元。

Scope of delivery 2.136.0305

| Qt. | Order no. | Description |
|-----|-----------|-------------|
|-----|-----------|-------------|

|       |            |                                                      |
|-------|------------|------------------------------------------------------|
| 1 PCS | 1.136.0305 | 用于 864 Robotic Balance Sample Processor 的离子交换器， 115V |
|-------|------------|------------------------------------------------------|

用于 864 Robotic Balance Sample Processor 的离子交换器（HAUG），带 115V 供电单元。



1 PCS


6.9013.000

Ionisation rod for ionizer

The ionisation rod is mounted with the holder (6.2057.160) to the 864 Robotic Balance Sample Processor.



# Optional accessories

| Order no.  | Description                                     |                                                                                     |
|------------|-------------------------------------------------|-------------------------------------------------------------------------------------|
| 6.2057.160 | 用于 864 Robotic Balance Sample Processor 的电离棒的支架 |  |

---