



826 pH mobile, 帶便攜箱

2.826.0110

便携式 pH 计，带 IrDA 接口。带便携箱和用于校正的多种附件，包括 Primatrode 6.0228.020。

Scope of delivery 2.826.0110

Qt.	Order no.	Description
1 PCS	1.826.0010	826 pH mobile

便携式 pH 计，带 IrDA 接口，包括背带和电池。



1 PCS 6.2050.000 pH计的背带
用于704, 826型pH计的背带。



1 PCS 6.2716.040 826 pH计的机盒
826 pH计及配件的机盒



1 PCS 6.1236.050 套管, SGJ 14/12 毫米
套管, SGJ 14/12 毫米, 聚乙烯



1 PCS 6.1446.000 SGJ 塞子 / B-14/ (15)



1 PCS

6.1613.020

试剂瓶 / 25毫升 / 带pH 7标签 (空瓶)



1 PCS

6.1613.030

试剂瓶 / 25毫升 / 带pH 4标签 (空瓶)



1 PCS

6.1614.000

洗瓶 / 250毫升



3 PCS

6.1614.030

塑料瓶 / 50毫升

塑料瓶



1 PCS 6.2133.000 5号碱性电池，1.5V，4只

每套包含4节1.5伏的5号碱性电池。



1 PCS 6.2307.230 缓冲溶液，pH 4，7与9

混合缓冲液，pH 4.00/7.00/9.00（25℃），装于单次使用的小袋中，
无色，盒装，每盒3 x 10 x 30 毫升



1 PCS 6.2308.050 电解液，3 mol/L KCl（50 mL）

电解液c (KCl) = 3 mol/L，50 mL（用于Ag/AgCl参比系统）



1 PCS 6.2621.140 六角扳手（Hexagon key），2.5毫米



1 PCS

6.2717.000

PP 杯, 100 mL

PP 聚丙烯制成的量杯, 100 mL。



Optional accessories

Order no.	Description
6.2325.000	pHit kit pH电极的维护套件该套件包括:- 50 mL 清洗液- 50 mL 3M KCl 溶液- 50 mL 储藏液- 2 储存容器- 使用说明

