

pH Lab Irda 230V EU

2.827.0014

Scope of delivery 2.827.0014

Qt.	Order no.	Description	
1 PCS	6.2013.010	定位环 用于直径为 10 mm 的支架杆。	
1 PCS	6.2016.050	支撑杆 / 300毫米	
1 PCS	6.2021.020	电极架 电极架，可用于4只电极与2只滴定头	

1 PCS

6.2161.010

插入式电源单元, 欧盟 230V/6V DC

电源, 用于827型 pH lab



1 PCS

6.2621.070

内六角扳手, 5 mm, 用于 IC Sample Processor

