

pH Lab Irda 230V UK

2.827.0019

Scope of delivery 2.827.0019

Qt.	Order no.	Description
1 PCS	6.2013.010	定位环 用于直径为 10 mm 的支架杆。
		
1 PCS	6.2016.050	支撑杆 / 300毫米
		
1 PCS	6.2021.020	电极架 电极架，可用于4只电极与2只滴定头
		

1 PCS

6.2621.070

内六角扳手，5 mm，用于 IC Sample Processor



1 PCS

6.2161.010

插入式电源单元，欧盟 230V/6V DC

电源，用于827型 pH lab

