



827 pH lab; 230 V, UK 帶 Primatrode

2.827.0119

实验室 pH 计，帶 IrDa 接口和 Primatrode 6.0228.010

Scope of delivery 2.827.0119

Qt.	Order no.	Description
1 PCS	6.2161.030	插入式电源单元，英国 230V/6V DC 电源，用于827型 pH lab



1 PCS

6.2013.010

定位环

用于直径为 10 mm 的支架杆。



1 PCS

6.2016.050

支撑杆 / 300毫米



1 PCS

6.2021.020

电极架

电极架，可用于4只电极与2只滴定头



1 PCS

6.2621.070

内六角扳手，5 mm，用于 IC Sample Processor



Optional accessories

Order no.	Description
-----------	-------------

6.2307.230 缓冲溶液，pH 4，7与9

混合缓冲液， pH 4.00/7.00/9.00（25℃）， 装于单次使用的小袋中，
无色， 盒装， 每盒3 x 10 x 30 毫升



6.2325.000 pHit kit

pH电极的维护套件该套件包括:- 50 mL 清洗液- 50 mL 3M KCl 溶液- 50
mL 储藏液- 2 储存容器- 使用说明

