



891 Professional Analog Out

2.891.0010

The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.

Scope of delivery 2.891.0010

Qt.	Order no.	Description
1 PCS	1.891.0010	891 Professional Analog Out

The 891 Professional Analog Out is an additional module for intelligent IC instruments. It allows to transfer analog signals of different chromatographic data to external data systems.



1 PCS

6.2128.200

电缆 891 模拟

用于连接到外部 AD 转换器上的模拟电缆



1 PCS

6.2156.070

谱峰思维系列离子色谱仪与扩展模块之间的连接电
缆，1 m

谱峰思维系列离子色谱仪与扩展模块之间的连接电缆

