

## Zubehör-Kit: OMNIS IC AS

6.05000.050


Scope of delivery 6.05000.050

| Order Number | Description               | Quantity |
|--------------|---------------------------|----------|
| 6.01816.010  | 硅胶软管, 8/12 mm, 2 m<br>排出管 | 2 PCS    |




---

|             |         |       |
|-------------|---------|-------|
| 6.02040.000 | AS Tray | 1 PCS |
|-------------|---------|-------|



---


|             |                            |       |
|-------------|----------------------------|-------|
| 6.02041.510 | OMNIS Probenrack 40 x 5 mL | 2 PCS |
|-------------|----------------------------|-------|



6.02061.120

**Tropfensammler, Extension**


1 PCS



6.02743.027

**Vial PS, 5 mL, 200 Stk.**

1 PCS



6.02743.030

**PP 瓶, 50 mL, 12 只**

1 PCS



50 mL 瓶由聚丙烯 (PP) 塑料制成, 可用于样品、标准、校验标准和冲洗液。在无菌和洁净室条件下进行生产和包装。这保证了较小的污染度。随机检查样品容器是否可渗出阴离子和阳离子。

6.02743.087

**Kappe PP, AD 11.5 mm, 200 Stk.**

1 PCS

