

Zubehör-Kit: OMNIS IC AS COOL


6.05000.060

Scope of delivery 6.05000.060


Order Number	Description	Quantity
6.01816.010	硅胶软管, 8/12 mm, 2 m 排出管	2 PCS



6.02041.510	OMNIS Probenrack 40 x 5 mL	2 PCS
-------------	----------------------------	-------




6.02061.120	Tropfensammler, Extension	1 PCS
-------------	---------------------------	-------



6.02743.027

Vial PS, 5 mL, 200 Stk.

1 PCS



6.02743.030

PP 瓶, 50 mL, 12 只

1 PCS




50 mL 瓶由聚丙烯 (PP) 塑料制成, 可用于样品、标准、校验标准和冲洗液。在无菌和洁净室条件下进行生产和包装。这保证了较小的污染度。随机检查样品容器是否可渗出阴离子和阳离子。

6.02743.087

Kappe PP, AD 11.5 mm, 200 Stk.

1 PCS



6.02040.010

OMNIS IC AS COOL Tray

1 PCS

