


OMNIS IC-Ausr. D-INJ, PP

6.05331.020

Scope of delivery 6.05331.020

Order Number	Description	Quantity
6.01826.020	Pumpschlauch, ID 1.0 mm, weiss - weiss	1 PCS



6.1831.180	PEEK Capillary 0.5 mm i.d., 3 m	1 PCS
------------	---------------------------------	-------



6.2744.034	2 个 Olive/UNF 10/32 两通	1 PCS
------------	------------------------	-------



连接压力螺丝和泵管。2 个。用于带有蠕动泵的离子色谱仪。

6.2744.070

压力螺丝（短）

1 PCS

较短的版本。带UNF 10/32接口。5件装。PEEK泵头的
备用件。

