



滴定头 / 管长0.4米

6.1537.010

Scope of delivery 6.1537.010

Qt.	Order no.	Description
1 PCS	6.1446.030	球塞 (Ball stopper) / B-14/ (15) 用于将滴定头可移动式地固定在滴定杯中



1 PCS

6.1543.060

滴定管头 / 螺纹 M6

滴定管头，螺纹 M6。与毛细管 6.1805.160 和可移动球塞 6.1446.030 一同组成整套的计量管头 6.1537.010。



1 PCS

6.1805.100

FEP管 / M6 / 40厘米

带光照保护和扭折保护

